

# In situ hybridization probe design and synthesis

## Reagents

### Nuclease Free H2O

HyClone HyPure Molecular Biology Grade Water. GE SH30538.02

### DNeasy Blood and Tissue Kit

Qiagen catalog #69504

### PCR Master Mix

Phusion Flash High-Fidelity PCR Master Mix. ThermoFisher F548S

### T7 RNA polymerase and buffer

ThermoFisher EP0112

### DIG-labelled NTP mix

Roche / Sigma 11277073910

### RNAse inhibitor

SUPERase-in, ThermoFisher AM 2694

### Linear polyacrylamide

VWR K548-1mL

### Formamide

Fisher, BP 227-500

## Design probes

We design 200-300bp probes, although we note that other labs design longer probes for fluorescent *in situ* hybridization experiments in *Drosophila* (<http://fly-fish.cabr.utoronto.ca/>).

For a given gene of interest, Identify a large exon using the Get Decorated FASTA tool from [flybase.org](http://flybase.org). Ensure that exon is present in all annotated transcripts using the UCSC Genome Browser.

Use Primer3Plus to select primers that will amplify a 200-300 bp region of the selected exon. Use the In Silico PCR tool from the UCSC Genome Browser to ensure that the selected pair will amplify a single region from the genome. If desired, use the In Silico PCR tool to screen suggested primer pairs from Primer3Plus for their ability to amplify the same region of DNA in multiple *Drosophila* species.

Order your selected primers using IDTDNA or another service, ensuring that the reverse primer (the primer that forms the 5' end of the antisense strand) has a T7 RNA polymerase binding sequence (5' TAATACGACTCACTATAG 3') on its 5' end.

## Amplify probe template

Isolate genomic DNA from your experimental strain using the (Qiagen) DNeasy kit.

PCR the probe template from genomic DNA. We use Phusion Flash Master Mix (Thermo Fisher) and associated protocol.

Run a small amount of the reaction (2uL + 8uL water) on an agarose gel to ensure that no off-target amplification has occurred. You should see a single band of the correct size.

## Synthesize probes

Use *in vitro* transcription to synthesize antisense RNA probes labelled with digoxigenin (DIG). We find that raw PCR products work well as template DNA.

	20 uL Reaction
<b>5X buffer</b>	4 uL
<b>10x NTP mix</b>	2 uL
<b>Template DNA (PCR reaction)</b>	2 uL
<b>RNAse Inhibitor</b>	1 uL
<b>T7 RNA polymerase</b>	1.5 uL
<b>Sigma water to final volume</b>	9.5 uL

Incubate at 37°C overnight. *We find that overnight incubation can improve yields over standard incubation times of 2-4 hours.*

If desired, run an agarose gel after incubation but prior to precipitation that compares equal volumes of transcription reaction and a 1:10 dilution of the template DNA. *This step checks the efficacy of the in vitro transcription reaction.*

### **Precipitate probes**

Add 1:10 volume (2uL) 3M NaOAc pH 5.2.

Add 2 volumes (44uL) 100% ethanol.

Add 2 uL linear polyacrylamide. *This decreases the amount of incubation time needed for ethanol precipitation.*

Incubate at -20°C for at least 30 minutes.

Spin at 4°C at max speed for 10 minutes. *A white pellet should be visible after this step*

Remove supernatant, and replace with 100uL 70% ethanol.

Spin at 4°C at max speed for 10 minutes.

If desired, repeat wash with 100% ethanol.

Remove supernatant, and air dry for 10 minutes.

Resuspend in 27 uL RNase-free water. Use 2 uL to quantify concentration and purity on a Nanodrop. Add 1 volume (25uL) formamide and store at -20°C.